

An A4 paper-sized blood vessel model that is easy to carry
 You can practice cerebrovascular access, thrombus retrieval, and aneurysm coil embolization.
 *Sold to medical institutions

Cerebrovascular smart model 3in1 (Catheter treatment training kit)

Joint development

Contracting
 experimental work


Products

Evaluation & testing

Academic networking


You can learn the following:

- Pass a guiding catheter from the femoral artery to the cerebral vessels
- Use a Simmons-shaped catheter from the brachial artery
- Simulated thrombus retrieval: Occlusion of the terminal part of the internal carotid artery and distal occlusion of the middle cerebral artery
- Catheter guidance and coil insertion into aneurysm




Design registration No. 1724174

Approach areas



Bovine arch



Aneurysm just past
 the bifurcation



Ledge at the origin
 of ophthalmic artery